## ActiMaris® – Wundspüllösung (300ml)

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## ActiMaris® – Wound irrigation solution (300ml)

### 1. Preliminary remark:

ActiMaris® wound irrigation solution provides outstanding cleaning on the tissue surface, both on biofilms and in hard-to-reach areas of the body.

### 2. Product profile and indications:

#### 2.1 For cleaning and irrigating, moistening and decontamination of

- a) acute wounds:
- mechanical wounds: skin lacerations, bites, cuts and abrasions and tears and crush wounds
  - postoperative wounds
- b) chronic wounds:
- decubitus ulcers
  - arterio-venous ulcers
  - diabetic ulcers
- c) thermal and chemical wounds:
- first- and second-degree burns
- d) entry portals of urological catheters, PEG tubes and drains

#### 2.2 Intraoperative cleaning and irrigation of wounds

e.g. in severe injuries of the extremities.

#### 2.3 For moistening bandages and wound dressings

such as compresses, gauze dressings, pads, sponges, gels, hydrofibres, alginates, hydrocolloids etc.; during dressing changes to loosen encrusted bandages or other encrusted wound dressings.

### 3. Use in general:

The wounds should always be irrigated and cleaned first with ActiMaris® wound irrigation solution so that easily detached coating is removed. Bandages or wound dressings soaked with ActiMaris® wound irrigation solution can be used. Because of the relatively low moisture binding, the dressing must be changed at least daily, and ideally several times a day. It should be used often until all coating and necrotic tissue can be removed easily and the wound is clean visually.

#### 4. Decontaminating cleaning of the skin, e.g. with a positive MRSA swab:

If necessary, decontamination of the entire part of the body or even the patient's entire body by washing with suitable means may be advisable.

#### 5. Use with encrusted dressings or when dressing changes are difficult:

Where dressings are difficult to remove, moistening the dressings with ActiMaris® wound irrigation solution is recommended until it is possible to separate the dressing gently without traumatising the wound surface.

#### 6. Tissue tolerability and biocompatibility:

Not cytotoxic; not irritating; not sensitising; low-pain, no inhibition of granulation and epithelialisation; assessed as safe dermatologically.

#### 7. Side effects:

None known at present. It cannot be ruled out that coloured textiles may be bleached after contact with ActiMaris®.

#### 8. Contraindications:

Since allergies can never be ruled out, ActiMaris® wound irrigation solution should not be used in the case of known allergies or suspected allergy to any of the ingredients.

#### 9. Restrictions on use:

Pregnancy and lactation: there is no evidence that the ingredients have mutagenic or embryotoxic effects. However, because of the lack of studies and long-term experience with pregnant and lactating women, ActiMaris® wound irrigation solution should be used in these cases only after careful medical assessment of the benefits and risks.

**Infants and toddlers:** because of inadequate long-term experience, ActiMaris® wound irrigation solution should be used in these cases only if highly indicated.

#### 10. General safety instructions:

Only for external use for skin wounds. Do not use for infusion or injection. Not to be taken orally. Use only intact and undamaged containers. Protect bottles from direct sunlight. Medical product. Keep out of reach of children.

#### 11. Composition / technical details:

Purified water, sea salt, activated oxygen

**Shelf life:** according to the expiry date at room temperature.

**Period of use:** 3 months after opening the bottle. Ensure that the container is always closed securely.

-  Lot number
-  Order number
-  Follow instructions for use
-  Use by: year and month
-  Manufacturer

**Date of information:** 01/2010